



# Hitachi 3D TEM

# Advantages

Using only **one TEM sample**, you are able to:

- Perform atomic-resolution imaging and electron diffraction without orientation limit.
- Quick 3D panoramic view for shape, dimension, location, and failure analysis of nanomaterials, nanocatalysts, thin film devices, and biological materials.
- Collect image data in the 0-180° angle range for electron tomography. No artifact for 3D reconstruction because of no missing data.
- 3D EDX and EELS chemical mapping.
- Prepare TEM samples using FIB and move to TEM and STEM instruments for characterization without having to relocate the tiny samples.

# System Requirements

## TEM: Hitachi transmission electron microscopes


- H-9500 300 kV high-resolution TEM for materials and biological sciences.
- H-7650 120 kV automated TEM for biological and beam-sensitive materials.
- HD-2300 200 kV dedicated STEM for semiconducting industry, materials science and biological science.

\*Please find more information at [www.hitachi-hita.com](http://www.hitachi-hita.com)

## Specimen Holder: Hitachi-patented 3D holder


- 0-360° rotation and  $\pm 15^\circ$  tilt under electron beam**
- Compatible with Hitachi Focus Ion Beam (FIB) instrument

## Hitachi-Patented Sample Holder: 360° Panoramic View



Movie available upon request

# Hitachi TEM: H-9500



---

# Special Functions for H-9500 TEM

**Option 1:** In-situ gas injection-heating, atomic resolution


**Option 2:** 360°-view structural and chemical analysis

# Hitachi TEM: H-7650

Operation  
Panel

Camera  
Chamber

Digital CCD  
Camera



W or LaB<sub>6</sub>  
Electron  
Source

5-Axis Stage

19" LCD  
Monitor

Control PC

---

# Special Features for H-7650 TEM


Automatic electron tomography

Automatic particle searching

**360°-view structural and chemical analysis**



# Hitachi STEM: HD-2300




# Special Features for HD-2300 STEM

- High through-put imaging + analysis tool
- SEM-like easy operation
- Complimentary BF, DF, SE Images
- EELS / EDS / Spectrum Imaging
- **3D structural and chemical characterization (Tomography reconstructions using Hitachi 3D holder)**
- Low-loss and Core-Edge fine structures
- Diffraction-based measurements

**\*Cs Correction STEM (HD-2700) is now available**

# Full Tilting-Angle-Range Electron Diffraction Using ONE TEM Sample



**H-9500 TEM**

T. Yaguchi et al, Japan. First prize at the photo contest of Japan microscopy society in 2005


**HITACHI**  
Inspire the Next

Hitachi High Technologies


# High-resolution Imaging Using Hitachi 3D Sample Holder

One sample for all-direction imaging and electron diffraction

H-9500  
300 kV TEM




High resolution TEM images (a)-(c) and electron diffraction patterns taken from three different directions of the same Si sample.





## 3D View of Nanomaterials Using Hitachi 3D Sample Holder



Instrument : Hitachi 200 kV dedicated STEM, HD-2300  
Sample rotation: 0 to 360°


[Movie available upon request](#)

# 3D Chemical Mapping Using Hitachi 3D Sample Holder

Mapping

## Ni-base Superalloy

Ni



Cr



Ti



100nm

Instrument : Hitachi 200 kV dedicated STEM, HD-2300

Probe size : 1 nm,

Probe current : 1 nA


[Movie available upon request](#)

# 3D Chemical Mapping Using Hitachi 3D Sample Holder

**Thin film device**

**Instrument: Hitachi 200 kV  
dedicated STEM, HD-2300**

Movie available  
upon request



---

## Comments and detailed information? Please contact:

Xiao Feng Zhang  
Hitachi High Technologies America, Inc.  
Electron Microscope Division  
5100 Franklin Drive  
Pleasanton, CA 94588  
USA

Tel: (925)-218-2814  
Fax: (925)-218-3230  
Email: [xiao.zhang@hitachi-hita.com](mailto:xiao.zhang@hitachi-hita.com)



[www.hitachi-hita.com](http://www.hitachi-hita.com)