


# S-3200N Scanning Microscope

## Introduction Manual



REV. 0

1/27/98

S3200INT-0

## S-3200N STRONG POINTS

### 1. COMPREHENSIVE VISUAL OPERATION

- Graphic color menu driven operation with mouse
- Graphic input and measurement
- Full keyboard operation (IBM Keyboard)

### 2. EXPANSIVE FRAME MEMORY

- 2048×1920 Pixel memory
- Electronic digital zoom

### 3. FULL AUTOMATED OPERATION

- Auto Start
- Auto Filament Saturation
- Auto Gun Alignment
- Auto Stigmation / Focus
- Auto Brightness & Contrast
- Auto Photo Mode
- Operating Condition Memory

### 4. NEW OBJECTIVE LENS DESIGN

- Resolution improvement
- Improvement of low kV image


### 5. X-ray MICROANALYSIS (NEW CHAMBER DESIGN)

- High X-ray take-off angle of 35°
- The specimen held at a short working distance of 15mm

### 6. VARIABLE PRESSURE VACUUM SYSTEM

### 7. VARIOUS SPECIMEN STAGE


- 100×50 large eucentric stage
- Cooling stage




## Installation Wire Connections

- Fore Pumps (2) connect to the column console in back -- CN100
- Air Compressor connects to the column console in back -- CN100
- Auto-transformer connects to the column console in back -- CN100
- Display power connects to the column console in back --CN100

CN100




- SE Detector (CN26) on the column connects to the Analog PCB located in the Display console
- Robinson Detector (CN300) on the column connects to the Analog PCB located in the Display console
- CN76 from Display console connects to the HVC2 PCB in the column console (HV Tank).
- CN36A from the Display console connects to the Evac PCB in the column console.
- CN112, CN114, & CN120 from the Display console connects to the Evac Seq.2 PCB in the column console.
- CN50 and CN51 from the Display console connects to the CDTB ( column distribution ) PCB in the column console.



New Design Units of S-3200N

Dimension of S-3200N (unit:mm)




# S-3200N.

## Multifunction Unit

### Multifunction Unit



#### Multifunction Mode Uses:

- ◆ Filament
- ◆ Gun Bias
- ◆ Raster Rotation
- ◆ Condenser Lens
- ◆ Line Position
- ◆ Spot Position
- ◆ Dynamic Focus
- ◆ Tilt Compensation





# S-3200N

## Signal Panels on Back of Display



# A F S   S y s t e m

Auto Filament Saturation System





C o n c e p t   o f   A F S   S y s t e m


# A G A ( Auto Gun Alignment )

Adjust gun alignment automatically so that the electron beam centers on the electron optics.

AGA MENU



# S-3200N Beam Current




Gun Bias=0

Filament Current=80 as read on the software icon.

.3-5KV=470K ohm Bias Res., 6-30KV=4M ohm Bias Res.

# S-3200N

## Specimen Vacuum vs Gun Vacuum




Ultimate Vacuum after 15 Hrs. =  $1.3 \times 10^{-6}$  Torr


Specimen Pump Down Time = 3min. 10sec.

Time to go to 1 PA from 270 PA = 1min. 30sec.

S - 3 2 0 0 N




Calculated Resolution of the S-3200N



Calculated Resolution of the S-3200N

Realization  
of Graphical  
User Interface


Clear Image  
by Memory Photo



Note (GDP :Graphic Data Processor  
DMAC:Direct Memory Access Controller  
LUT :Look-Up Table)

Configuration of Image Memory Unit

Chamber of S-3200N (6 inch)



Name	Model	Part No.	Mounting Port
Specimen tensil stage:200kg	S-5757	49E-0142	Ⓐ
Z variable cryogenic system		49E-0232	Ⓐ
Specimen manipulator	S-5744	50E-0240	Ⓒ
Cathodoluminescence detector		48E-0207	Ⓓ or Ⓛ
Quartered-semiconductor backscattered electron detector		50E-0100	Ⓒ or Ⓛ
Farady cup device	S-5766	48E-0249	Ⓒ or Ⓚ
Energy dispersive x-ray detector			Ⓑ

## S P E C I F I C A T I O N S

1. Resolution	3.5nm(at high vacuum mode) 5.5nm(at low vacuum mode)
2. Magnification	20 to 300,000×
3. Electron Optics	(1) Accelerating voltage 0.3 to 10kV(0.1kV step) 10 to 30kV(1kV step) (2) Electron gun Pre-centered hairpin type filament (3) Bias Self-bias exchangeable + continuously variable (4) Gun alignment 2-stage electromagnetic alignment (5) Lens system 3-lens reduction type(super-conical objective lens) (6) Objective lens aperture Movable aperture (4 openings selectable and finely adjustable) (7) Stigmator coil 8-pole electromagnetic type (8) Deflection coil 2-stage electromagnetic deflection type(X,Y)
4. Specimen Size	150mm dia. (MAX)
5. Specimen Stage	Any one of the specimen stages(1),(2),(3) and (4) below is equipped. (1) Super eucentric stage X-traverse: 0 to 32mm (continuous) Y-traverse: 0 to 32mm (continuous) Tilt angle: -90° to 90° (continuous) Rotation angle: 360° (continuous) Z-traverse: 5 to 35mm (continuous) 10 to 35mm (with BSE detector) (2) Large-size standard stage X-traverse: 0 to 80mm (continuous) Y-traverse: 0 to 40mm (continuous) Tilt angle: -20° to 90° (continuous) Rotation angle: 360° (continuous) Z-traverse: 5 to 35mm (continuous) 10 to 35mm (with BSE detector) (3) Large-size eucentric stage X-traverse: 0 to 100mm (continuous) Y-traverse: 0 to 50mm (continuous) Tilt angle: 0° to 60° (continuous) Rotation angle: 360° (continuous) Z-traverse: 5 to 35mm (continuous) 10 to 35mm (with BSE detector)

**9. Evacuation System**

- (1) Sequence control  
Full-auto pneumatic valve system
- (2) Low vacuum control  
Automatic (real-time vacuum feed-back system)
- (3) Vacuum detection  
With Pirani gauge
- (4) Ultimate vacuum  
 $6 \times 10^{-4}$  Pa
- (5) Low vacuum setting & indication  
On menu screen
- (6) Low vacuum range  
1 to 270 Pa
- (7) Evacuation time  
About 3 min.
- (8) Vacuum pump  
Oil rotary pump (160 l/min) × 2  
Oil diffusion pump (570 l/sec) × 1
- (9) Air compressor  
Take Off Angle: 35° (at WD15 mm)

**10. X-ray Analysis**

Position

**11. Video Out**

**12. Safety Device**

**13. Water Supply  
and Drain**

**14. Dimensions**

NTSC-TV signal (B/W)

Protection against power, water and vacuum failures

(1) Water flow rate

1 to 1.5 l/min

(2) Water pressure

50 to 100 kPa (0.5 to 1 kgf/cm<sup>2</sup>)

(3) Water temperature

10 to 20°C

(1) Main console

700 (W) × 700 (D) × 1500 (H) mm, 320 kg

(2) Display unit

1000 (W) × 900 (D) × 1104 (H) mm, 180 kg

(4) Cool stage  
 X-traverse: 0 to 15mm (continuous)  
 Y-traverse: 0 to 15mm (continuous)  
 Tilt angle: -45° to 45° (continuous)  
 Rotation angle: unavailable  
 Z-traverse: 10 to 35mm  
 Temperature range: -10 to 10°C

6. Detector  
 (1) Secondary electron detector  
 (2) High-sensitivity (Robinson type) backscattered electron detector

7. Operation  
 Graphical User Interface with mouse  
 (operation with rotary knob and keyboard available)

8. Image Display  
 (1) Principle  
 Flicker-free image display with built-in image memory  
 (2) Frame memory  
 $512 \times 480 \text{ pixel} \times 16 \text{ bit}$   
 $2048 \times 1920 \text{ pixel} \times 8 \text{ bit}$   
 (3) Viewing CRT  
 14-inch TV monitor (Color)  
 (4) Photo CRT  
 7-inch Ultrahigh resolution CRT (B/W)  
 (5) Scanning mode  
 TV scan  
 SLOW scan  
 PHOTO scan  
 Selected area scan  
 Waveform monitor/signal monitor  
 TWIN scan  
 Dual-Mag scan  
 (6) Scan speed  
 TV:  $\text{Y}=0.033 \text{ sec } (\text{X}=0.064 \text{ msec}, \text{LINES}=500)$   
 SLOW1:  $\text{Y}=0.35 \text{ sec } (\text{X}=0.7 \text{ msec}, \text{LINES}=500)$   
 SLOW2:  $\text{Y}=2 \text{ sec } (\text{X}=4 \text{ msec}, \text{LINES}=500)$   
 SLOW3:  $\text{Y}=10/8.4 \text{ sec } (\text{X}=20/16.7 \text{ msec}, \text{LINES}=500)$  \*  
 SLOW4:  $\text{Y}=20/25 \text{ sec } (\text{X}=40/50 \text{ msec}, \text{LINES}=500)$  \*  
 S.A:  $\text{Y}=40 \text{ msec } (\text{X}=0.35 \text{ msec}, \text{LINES}=125)$   
 PHOTO1:  $\text{Y}=40/34 \text{ sec } (\text{X}=20/16.7 \text{ msec}, \text{LINES}=2000)$  \*  
 PHOTO2:  $\text{Y}=80/100 \text{ sec } (\text{X}=40/50 \text{ msec}, \text{LINES}=2000)$  \*  
 PHOTO3:  $\text{Y}=200/200 \text{ sec } (\text{X}=100/100 \text{ msec}, \text{LINES}=2000)$  \*  
 PHOTO4:  $\text{Y}=400/400 \text{ sec } (\text{X}=200/200 \text{ msec}, \text{LINES}=2000)$  \*  
 \*Power-synchronous scan: 50/60Hz region

(4) Cool stage  
 X-traverse: 0 to 15mm (continuous)  
 Y-traverse: 0 to 15mm (continuous)  
 Tilt angle: -45° to 45° (continuous)  
 Rotation angle: unavailable  
 Z-traverse: 10 to 35mm  
 Temperature range: -10 to 10°C

**6. Detector**  
 (1) Secondary electron detector  
 (2) High-sensitivity (Robinson type) backscattered electron detector

**7. Operation**  
 Graphical User Interface with mouse  
 (operation with rotary knob and keyboard available)

**8. Image Display**  
 (1) Principle  
 Flicker-free image display with built-in image memory  
 (2) Frame memory  
 512 × 480 pixel × 16bit  
 2048 × 1920 pixel × 8bit  
 (3) Viewing CRT  
 14-inch TV monitor (Color)  
 (4) Photo CRT  
 7-inch Ultrahigh resolution CRT (B/W)  
 (5) Scanning mode  
 TV scan  
 SLOW scan  
 PHOTO scan  
 Selected area scan  
 Waveform monitor/signal monitor  
 TWIN scan  
 Dual-Mag scan  
 (6) Scan speed  
 TV: Y=0.033sec (X=0.064 msec, LINES=500)  
 SLOW1: Y=0.35sec (X=0.7 msec, LINES=500)  
 SLOW2: Y=2sec (X=4 msec, LINES=500)  
 SLOW3: Y=10/8.4sec (X=20/16.7 msec, LINES=500) \*  
 SLOW4: Y=20/25sec (X=40/50 msec, LINES=500) \*  
 S.A: Y=40 msec (X=0.35 msec, LINES=125)  
 PHOTO1: Y=40/34sec (X=20/16.7 msec, LINES=2000) \*  
 PHOTO2: Y=80/100sec (X=40/50 msec, LINES=2000) \*  
 PHOTO3: Y=200/200sec (X=100/100 msec, LINES=2000) \*  
 PHOTO4: Y=400/400sec (X=200/200 msec, LINES=2000) \*  
 \*Power-synchronous scan: 50/60Hz region

(7) Auto mode

Auto filament saturation(Auto emmision current)  
Auto gun alignment  
Auto focus  
Auto stigmator  
Auto brightness/contrast control  
Auto start(HV-on→ABC→AFC)  
Auto photo(1:AFC→ABC→PHOTO)  
(2:ABC→PHOTO)

(8) Signal processing

Recursive filter(TV,SLOW1 and SLOW2 scans)  
Image integration(TV,SLOW1 and SLOW2 scans)  
Brightness conversion  
Gamma control(digital)  
Polarity reversion(digital)  
Histogram  
Quartered frame copying  
Digital zoom

(9) Data entry

Arbitrary alphanumeric characters and symbols  
overlayable at any image position through keyboard

(10) Auto data display

Accelerating voltage,magnification,micron scale,  
micron value,film number,WD value,time,date and  
low vacuum set value printable on film

(11) Operation assist functions

Preset magnification  
Wobbler mode  
Registration of column conditions  
(15 conditions×50 sets)  
Focus search mode  
Focus memory  
Polarity reversion(analog)  
Differential image(analog)  
Gamma control(analog)  
Setting brightness/contrast for ABC  
Setting brightness/contrast for photo  
Setting film ASA  
Setting photo ratio  
Scan indicator  
Filament image  
Display of condition table

- 9. Evacuation System**
- (1) Sequence control  
Full-auto pneumatic valve system
  - (2) Low vacuum control  
Automatic(real-time vacuum feed-back system)
  - (3) Vacuum detection  
With Pirani gauge
  - (4) Ultimate vacuum  
 $6 \times 10^{-4}$ Pa
  - (5) Low vacuum setting & indication  
On menu screen
  - (6) Low vacuum range  
1 to 270Pa
  - (7) Evacuation time  
About 3min.
  - (8) Vacuum pump  
Oil rotary pump(160l/min) × 2  
Oil diffusion pump(570l/sec) × 1
  - (9) Air compressor  
Take Off Angle:35° (at WD15mm)
- 10. X-ray Analysis**
- Position
- 11. Video Out**
- NTSC-TV signal(B/W)
- 12. Safety Device**
- Protection against power, water and vacuum failures
- 13. Water Supply and Drain**
- (1) Water flow rate  
1 to 1.5l/min
  - (2) Water pressure  
50 to 100kPa(0.5 to 1kgf/cm<sup>2</sup>)
  - (3) Water temperature  
10 to 20°C
- 14. Dimensions**
- (1) Main console  
700(W) × 700(D) × 1500(H)mm, 320kg
  - (2) Display unit  
1000(W) × 900(D) × 1104(H)mm, 180kg

## OPTIONAL ACCESSORIES

### (1) Specimen Stages

Part No.	Article	Description
49E-0142	S-5757 200kg specimen tensile stage	Tensile load: 200kg max. Exclusive goniometer stage Order tensile load measuring device and stage adapter together with this stage. Recorder (or voltmeter) is to be prepared by customer.
47E-5224	Tensile load measuring device	For S-5757 Semiconductor load converter and DC amplifier
50E-3051	Cool stage	Stage + control unit -10~10°C (1°C step)
49E-0232	Z variable cryogenic system	For observing frozen specimen Charge for installation is included in price of this system. Oil rotary pump is to be prepared by customer.
47E-5229	S-4078 Evaporation unit	For cryogenic system Only gold coating possible.
50E-0410	Stage motor drive unit S II	For use with 5" stage. Order together with touch sensor (P/N 50E-0426)
50E-0411	Stage motor drive unit S II	For use with large-sized eucentric stage. Order together with touch sensor (P/N 50E-0426)
50E-0426	Touch sensor	For use with S I and S II
49E-0645	Stage motor drive unit S I	For use with 5" stage. Order together with touch sensor (P/N 50E-0426)
50E-0240	Specimen manipulator	
47E-0244	Specimen holder (four 15mm dia. specimens)	For 5" stage and large-sized eucentric stages
47E-0245	Specimen holder (twenty-one 6mm dia. specimens)	

(2) Display Units

Part No.	Article	Description
50E-0202	English/Japanese menu changeable function	Order together with main unit.
50E-0204	Raster rotation/dynamic focus/tilt compensation unit	
50E-0205	Image filing 3.5-inch magneto-optical disk device(type 230MB)	(1) 3.5-inch magneto-optical disk device (2) One 3.5-inch magneto-optical disk (3) Cleaning kit
50E-0206	3.5-inch magneto-optical disk(type 230MB)	For P/N50E-0205(preformatted)
50E-0207	RS-232C communication interface	Order together with main unit.
50E-0208	Simple measurement /graphic input function	Order together with main unit.
50E-0209	Pseudo color function	Order together with main unit.

(3) Detectors

Part No.	Article	Description	
48E-0204	S-6542 Absorbed electron detector	For 5"stage Preamp only	Each detector includes no video amp. So order P/N50E -0201 together (one unit required for up to 3 detectors).
48E-0326	24-pin IC holder	Usable for 5"stage only	
48E-0207	S-5744 Cathodeluminescence detector	Detector+Preamp Order P/N50E-0200 together.	
50E-0100	Quartered-semiconductor backscattered electron detector	For S-3200N only	
50E-0200	Photomultiplier power supply	For S-5744	
50E-0201	Video amp unit	For S-6542/S-5744/quartered-semiconductor backscattered electron detector	

(4) Camera and Film

Part No.	Article	Description
50E-0250	Camera body	S-5080A Camera and X1adapter
50E-0251	6×7 Roll film unit	Combined with camera body. Consists of roll film holder, X0.6adapter and focusing glass.
50E-0252	107(105) Polaroid unit	Combined with camera body. Consists of 107 Polaroid film holder, X0.8 adapter and focusing glass.
50E-0254	4"×5" Polaroid unit (8 exposures)	Combined with camera body. Consists of 4"×5" Polaroid film holder (8 exposures) and focusing glass.
50E-0253	4"×5" Polaroid unit (1 exposure)	Combined with camera body. Consists of 4"×5" Polaroid film holder (1 exposure) and focusing glass.

Unit Name		Camera Body	6×7 Roll	107 Polar.	4"×5" Polaroid
Camera Unit	P/N	50E-0250	50E-0251	50E-0252	50E-0254 or 50E-0253
Purchase of 6×7 only		○	○		
Purchase of 107 only		○		○	
Purchase of 4"×5" only		○			○
Purchase of 6×7&107		○	○	○	
Purchase of 6×7&4"×5"		○	○		○
Purchase of 107&4"×5"		○		○	○
Purchase of 6×7, 107 & 4"×5"		○	○	○	○

(5) Items for X-ray Elemental Analysis

Part No.	Article	Description
50E-0210	X-ray mode function	
48E-0249	Farady cup unit	
50E-0340	Movable objective aperture plate for EDX	Used for excluding system peak
48E-0446	Antivibration rubber pad for EDX	Necessary for mounting EDX analyzer.
50E-0212	DBC interface unit	


HS  
6.12.72  
富澤

主要基板 VR・SW・JP 設定表

NOT USE

基板名	SW・VR名	内容
SH-SBC (150E2110)	SW1 A	0 : N-SEM & W 1 : N-SEM & LaBe 2 : H-SEM & W 3 : H-SEM & LaBe 4 : UPD & W 5 : UPD & LaBe 6 : 7 :
	SW1 B	0 : 1 : 2 : 3 : 4 : 5 : 6 : 7 :
	SW2 JAPANESE Ver	OFF 1 → 日本語メニュー 2 3 → X線モード無し 4 → 疑似カラー無し 5 6 7 → 日本語/英語切替なし 8 → 日本語/英語切替あり
NO. MEASURE / Graphic	SW3	1 R. R無し(Nope) 2 RS-232C無し( ) 3 SCSI無し( ) 4 5 低真空0-270Pa 6 アニュラ(BSE2)無し(Nope) 7 MOD無し( ) 8 クラフィック無し
	SW4	1 2 3 4 5 6 7 8
	JP1	SHORT : WDT有り OPEN : WDT無し(デバッグ用)
	JP2	SHORT : プザー音大 OPEN : プザー音小
	ANALOG (150E2102)	VR1 : SLOW3,4:SG-X CYCLE(TX=4ms, TP3) VR2 : SLOW2 : SG-X CYCLE(TX=4ms, TP3) VR3 : SLOW1 : SG-X CYCLE(TX=0.7ms, TP3) VR4 : SLOW1-4:SG-Y AMPLITUDE(9.12V, TP4) VR5 : SPLIT・DUALMAG:SG-X AMPLITUDE(10.8V, TP9) VR6 : AFC-SCAN:X AMPLITUDE(12V, TP5) VR7 : X-RAY MODE:X POSITION(V=-10.9V, TP9) VR8 : AFC-SCAN:Y AMPLITUDE(8.8V, TP6) VR9 : X-RAY MODE:Y POSITION(V=-9.6V, TP10) VR10 : PHOTO-CRT:X AMPLITUDE VR11 : PHOTO-CRT:X POSITION VR12 : PHOTO-CRT:Y AMPLITUDE VR13 : PHOTO-CRT:Y POSITION VR14 : X-RAY MODE:X AMPLITUDE(V+=10V, TP9) VR15 : X-RAY MODE:Y AMPLITUDE(V+=8.7V, TP10)


(ANALOG)	VR 1 6 VR 1 7 VR 1 8 VR 1 9 VR 2 0 VR 2 1 VR 2 2 VR 2 3 VR 2 4 VR 2 5 VR 2 6 VR 2 7 VR 2 8 VR 2 9 VR 3 0 VR 3 1 VR 3 2 VR 3 3 VR 3 4 VR 4 0	(MAG-X:UPD) MAG-X ORTHOGONAL ADJ. (ORTHOGONAL ADJ.:UPD) (MAG-X:UPD) MAG-Y TV-SCAN:X AMPLITUDE TV-SCAN:X POSITION TV-SCAN:Y AMPLITUDE TV-SCAN:Y POSITION ABC-CONTRAST:ADC LEVEL ABC-BRIGHTNESS:ADC LEVEL LINE-ANALYSIS1 AMPLITUDE LINE-ANALYSIS1 POSITION LINE-ANALYSIS2 AMPLITUDE LINE-ANALYSIS2 POSITION TV-SCAN:MEMORY AREA PHOTO-CRT BRIGHTNESS(COURSE) PHOTO-CRT BRIGHTNESS(FINE) BSE1 SIGNAL:GAIN
	J P 1 J P 2 J P 3 J P 4 J P 5 J P 6 J P 7	SHOR T : DBC 無し OPEN : DBC 有り (X - SCAN) SHOR T : DBC 無し OPEN : DBC 有り (Y - SCAN) SHOR T : R. R 無し OPEN : R. R 有り (X - SCAN) SHOR T : R. R 無し OPEN : R. R 有り (Y - SCAN) SHOR T : DBC 無し OPEN : DBC 有り (BLK) Composite BC (Analog) SHOR T : DBC 無し OPEN : DBC 有り (VIDEO) ALWAYS SHOR T
DEF / LENS (150E2103)	VR 1 VR 3 VR 4 VR 5 VR 6 VR 7 VR 8 VR 10	OBJ - SPAN COND 1 - SPAN COND 2 - SPAN TV - X : MAG CENTER SLOW - X : MAG CENTER TV - Y : MAG CENTER SLOW - Y : MAG CENTER REFERENCE : +10V (TP 1)
HVC 2 (150E2125)	VR 1 VR 2 VR 3	ACC. VOLTAGE ADJ. FIL. CURRENT ADJ. FREQUENCY (32 μ s)
SEQ 2 (250E2132)	VR 1 VR 2 VR 3 J P 1 J P 2 J P 3 J P 4 J P 5 J P 6 J P 7 J P 8	GUN-VACUUM LEVEL ADJ. CHAMBER-VACUUM LEVEL ADJ. SPECIMEN-CHANGER VACUUM ADJ. 1 - 2 SHOR T : HIGH 2 - 3 SHOR T : NATURAL SHOR T : UTW 無し OPEN : UTW 有り 1 - 2 SHOR T : W 電子銃 2 - 3 SHOR T : LaB <sub>6</sub> 電子銃 1 - 2 SHOR T : W 電子銃 2 - 3 SHOR T : LaB <sub>6</sub> 電子銃 ALWAYS SHOR T SHOR T : アニュラー無し OPEN : アニュラー有り SHOR T : ROBINSON 無し OPEN : ROBINSON 有り SHOR T : 試料交換装置無し OPEN : 試料交換装置有り



14-DEC-94

000033 WD12mm 1.0kV x450 100μm

Accelerating voltage  
: 1 kV  
Specimen  
:




14-DEC-94


000030 WD12mm 1.0kV x1.5k 20μm

Accelerating voltage  
: 1 kV  
Specimen  
: Bronchus of mouse


Improvement of Low kV



Digital ZOOM (1)




ZOOM( x 4 )



High Resolution Image (2048Pixel)

Digital ZOOM ( 2 )




HITACHI

S - 3 2 0 0 N

# Analog PCB

Arrangement for Gain 1-0

SG MODE Section



$I\bar{C}\bar{E}7$   
 $E\bar{D}$   
 $E\bar{9}$

SG MODE Section

(1) Functional block diagram

Mod 0 = Slow X  
 1 = Split X  
 2 = Line Analysis & Slow (LA slow)  
 3 = X ADD - 10V (summing)  
 4 = Slow Y  
 5 = Y ADD - 10V  
 6 = SCAN TO PCRT.

$T\bar{X}$  = Selected Area Analysis (Translational SAA)  
 $M\bar{S}\bar{A}\bar{A}$  = Mode Selected AREA ANALYSIS.

T-Y TRIANGLE Selected Analysis.  
 C-X = Autofocus  
 CIRCLE SCAN

(2) SG section

① SLOW-X

a) Scan speed

- The power synchronization speed (20 to 200ms) is controlled by controlling the integrator input current via 8 bit DAC. For 0.72 and 4ms, apply a current by a separate circuit.
- For the reduce area scan, limit the amplitude to 1/2.


X interval (50/60Hz)	Running time	Dwell time	AC Synchronization
0.72ms	0.62ms	0.1ms	×
4ms	3.9ms	0.1ms	×
20/16.7ms	16/12.7ms	4ms	○
40/50.1ms	36/46.1ms	4ms	○
100ms	96ms	4ms	○
200ms	196ms	4ms	○
Reduce	0.35ms	0.1ms	×

b) Control function

- X RESET : Return to start point and stop
- X FREE : Free run
- ESYNC X : AC sync ON/OFF

$$\left. \begin{array}{l} \\ \\ \end{array} \right\} I \in H_3 : \text{Antcos}(1/4)$$


c) Output voltage (TP3)



② SLOW-Y

a) Raster line quantity control system

- Control the number of lines by counting X-BLK signal and resetting when a set number of lines is reached. Y Line Counter M<sub>5</sub>, M<sub>6</sub>, K<sub>5</sub> Analog



480L  
(2/4)


$$\text{Number of raster lines is: } N = \text{DATA} \times 16$$

(DATA : Y line data, DY7~0)


b) Amplitude control

- The integrating circuit integrates X-BLK pulses at a constant gain.

It is followed by an amplifier and its gain is changed according to the number of lines. Therefore, the output amplitude obtained is the same at any number of lines.



YLO = Photo 1, 2, 3, 4 & MEMORY = 19  
YLINE1 = (NOT USED) - 1000 Line Scan  
YLINE2 = 480 Lines Slow Scan 1, 2, 3, 4  
YLINE3 = ABC 9 AND 4 32 Lines



Binney Hex  
Z 16

DY70

DY6	1	7
5	,	
4	1	
<u>3</u>	-	1
<u>2</u>	-	0
<u>1</u>	-	0
0		8

No. of X lines	Y LINE data	Amplifier gain	Use	Output swing
1920 <del>2048</del>	78 [80]h	x1	Direct photo	+10v 9V
960 <del>1024</del>	3C [40]h	x2	Memory photo	+10v 9V
480 <del>512</del>	1E [20]h	x4	Slow	+10v 9V
256	[10]h	x4	Monitor (waveform)	5V
128	[08]h	x4	Reduce	2.5V
32	[02]h	x64	ABCC	10V


### c) Control function

- Y Reset: Return to start point and stop.
- Y FREE : Free run

Time Display - IMAGE only

### Analog 4/4 ③ TV-X/TV-Y - Column Scan

- a) Miller integrating type. Synchronized with TDISP(X), VD signal generated by the image memory (IMEM1).



### b) Control function


Set to ON/OFF by software program.

### c) Amplitude

Same as SLOW except that X over-scans beyond the range of 0 to 10v.  
(Adjust to match the SLOW image.)

### ④ C-X/C-Y

By sin/cos oscillator and integrator, outputs  $\sin t$  and  $\sin 2t$ .



At every interval, delivers a trigger output (CST). Used for starting AFC cycle.


### ⑤ T-X/T-Y


Generates a triangle wave of interval 4.8ms(X)/0.2s(Y) (for SAA).


Analog  
4/4

VIDEO AMP Section


(1) Functional block diagram







## PHOTO SIZE



### PHOTO

- PHOTO CRT DEF. AREA :  $126 \times 114\text{mm}$  ( $528 \times 480\text{pixel}$ )
- IMAGE MEMORY AREA :  $122 \times 114\text{mm}$  ( $512 \times 480\text{pixel}$ )
- CRT MONITOR VIEW AREA :  $122 \times 95\text{mm}$  ( $512 \times 400\text{pixel}$ )
- PHOTO FILM AREA :  $114 \times 89\text{mm}$  ( $479 \times 374\text{pixel}$ )
- ADJUSTMENT BOX :  $112 \times 87\text{mm}$  ( $470 \times 365\text{pixel}$ )

### VIEW MONITOR

- CRT MONITOR VIEW AREA :  $195 \times 152\text{mm}$  ( $512 \times 400\text{pixel}$ )
- CRT MONITOR MENU SIZE :  $244 \times 186\text{mm}$  ( $640 \times 480\text{pixel}$ )


VIEW RATIO of MONITOR :  $\times 1.6$

PIXEL ASPECT : 1 : 1

PHOTO PIXEL RESOLUTION : 4.2dot/mm (107DPI)


# Maintenance Menu 1

The MOUSE CURSOL move to left bottom edge on monitor, and push left button.




## Maintenance Menu 2

The menu display condition of mode switch on the SH-SBC.



( ) - NOT USED

PCB NAME	SW, VR NAME	FUNCTION
ANALOG	V R 1	SLOW3,4:SG-X CYCLE( $T_{hx}=4\text{ ms}$ ,TP3)
	V R 2	SLOW2 : SG-X CYCLE( $T_x=4\text{ ms}$ ,TP3)
	V R 3	SLOW1 : SG-X CYCLE( $T_x=0.7\text{ ms}$ ,TP3)
	V R 4	SLOW1-4:SG-Y AMPLITUDE(9.1V,TP4)
	V R 5	SPLIT·DUALMAG:SG-X AMPLITUDE
	V R 6	$\infty$ - SCAN:X AMPLITUDE(12V,TP5)
	V R 7	X-RAY MODE:X POSITION
	V R 8	$\infty$ - SCAN:Y AMPLITUDE(8.8V,TP6)
	V R 9	X-RAY MODE:Y POSITION
	V R 1 0	PHOTO-CRT:X WIDTH
	V R 1 1	PHOTO-CRT:X POSITION
	V R 1 2	PHOTO-CRT:Y WIDTH
	V R 1 3	PHOTO-CRT:Y POSITION
	V R 1 4	X-RAY MODE:X AMPLITUDE
	V R 1 5	X-RAY MODE:Y AMPLITUDE
	V R 1 6	(MAG-Y:UPD) - upper deflector $\leftarrow$ NOT USED
	V R 1 7	MAG-X(COLUMN SCAN X WIDTH) CW SMALLER
	V R 1 8	COLUMN SCAN RIGHT ANGLE
	V R 1 9	(COLUMN SCAN RIGHT ANGLE:UPD) -
	V R 2 0	(MAG-X:UPD)
	V R 2 1	MAG-Y(COLUMN SCAN Y WIDTH) CW SMALLER
	V R 2 2	TV-SCAN:X AMPLITUDE
	V R 2 3	TV-SCAN:X POSITION
	V R 2 4	TV-SCAN:Y AMPLITUDE
	V R 2 5	TV-SCAN:Y POSITION
	V R 2 6	ABC-CONTRAST:ADC LEVEL
	V R 2 7	ABC-BRIGHTNESS:ADC LEVEL
	V R 2 8	LINE-ANALYSIS1 AMPLITUDE
	V R 2 9	LINE-ANALYSIS1 POSITION
	V R 3 0	LINE-ANALYSIS2 AMPLITUDE
	V R 3 1	LINE-ANALYSIS2 POSITION
	V R 3 2	TV-SCAN:MEMORY AREA( $T_{hx}=10\text{ }\mu\text{s}$ ,TP24)
	V R 3 3	PHOTO-CRT BRIGHTNESS(COURSE)
	V R 3 4	PHOTO-CRT BRIGHTNESS(FINE)
	V R 4 0	BSE1 SIGNAL:GAIN
	J P 1	SHORT: NONE OPEN : DBC (X-SCAN)
	J P 2	SHORT: NONE OPEN : DBC (Y-SCAN)
	J P 3	SHORT: NONE OPEN : RASTER ROTATION (X-SCAN)
	J P 4	SHORT: NONE OPEN : RASTER ROTATION (Y-SCAN)
	J P 5	SHORT: NONE OPEN : DBC (BLK) composite
	J P 6	SHORT: NONE OPEN : DBC (VIDEO)
	J P 7	ALWAYS SHORT




NOTE • at DBC connected, remove JumperPin JP1, JP2, JP6, JP5  
• at Raster Rotation PCB fixed, remove JumperPin JP3, JP4.


Reserve

ANALOG PCB

PCB NAME	SW, VR NAME	FUNCTION
DEF / LENS	VR 1	OBJ-SPAN
	VR 3	COND1-SPAN
	VR 4	COND2-SPAN
	VR 5	TV-X : MAG CENTER — zoom & focus MAG ADJ.
	VR 6	SLOW-X : MAG CENTER
	VR 7	TV-Y : MAG CENTER — zoom to 350 MAG ADJ.
	VR 8	SLOW-Y : MAG CENTER
	VR 9	DYNAMIC FOCUS CENTER
	VR 10	REFERENCE : +10V(TP1)



DEF/LENS PCB





DCPS PCB

PCB NAME	SW, VR NAME	FUNCTION	
SH-SBC	SW 1 A NORMAL Position = 0	0 : N-SEM & W Tungsten 1 : N-SEM & LaB <sub>6</sub> 2 : H-SEM & W 3 : H-SEM & LaB <sub>6</sub> 4 : UPD & W 5 : UPD & LaB <sub>6</sub> 6 : 7 :	8 : 9 : A : B : C : D : E : F :
	SW 1 B NORMAL Position = 0	0 : 1 : 2 : 3 : 4 : 5 : 6 : 7 :	8 : 9 : A : B : C : D : E : F :
	SW 2 USA	OFF 1 JAPANESE 2 NTSC 3 NONE 4 NONE 5 RESERVE 6 RESERVE 7 8 NONE	ON ENGLISH PAL X-RAY MODE PSEUDO COLOR LANGUAGE CHANGE RR, DF, TC RS-232C SCSI 0-530Pa BSE2 ELEMENT MOD GRAPHIC
	SW 3 USA	1 NONE 2 NONE 3 NONE 4 RESERVE 5 0-270Pa 6 NONE 7 NONE 8 NONE	RR = RASTER Rotate DF = Dynamic focus TC = TILT Compens. 0-530Pa BSE2 ELEMENT MOD GRAPHIC
	SW 4 USA NORMAL - System 2 off	1 RESERVE 2 RESERVE 3 RESERVE 4 5 RESERVE 6 7 SYSTEM 1 OFF 8 SYSTEM 2 OFF	SYSTEM 1 ON SYSTEM 2 ON ONLY for SERVICE ADJUST.
	J P 1	SHORT : WDT NORMAL (AUTO RESET) OPEN : WDT NONE (for DEBUG)	
	J P 2	SHORT : BUZZER LOUD OPEN : BUZZER LOW	


System #1 = NO MEANING

System #2 = NO MEANING  
System 2 on - only for Service Engneese Adjust.

System #2 = ON MEANS SEBCAS,




PCB NAME	SW, VR NAME	FUNCTION
HVC 2	VR 1 VR 2 VR 3	ACC. VOLTAGE ADJUST FILAMENT CURRENT ADJUST FREQUENCY (32 $\mu$ S)





HVC 2 PCB

PCB NAME	SW, VR NAME	FUNCTION
SEQ 2	VR 1	GUN-VACUUM LEVEL ADJ.
	VR 2	CHAMBER-VACUUM LEVEL ADJ.
	VR 3	SPECIMEN-CHANGER VACUUM ADJ.
	J P 1	1-2 SHORT : HIGH 2-3 SHORT : NATURAL
	J P 2	SHORT: NONE OPEN : UTW
	J P 3	1-2 SHORT: W-GUN 2-3 HSORT: LaB6-GUN
	J P 4	1-2 SHORT: W-GUN 2-3 SHORT: LaB6-GUN
	J P 5	ALWAYS SHORT
	J P 6	SHORT: NONE OPEN : BSE2
	J P 7	SHORT: NONE OPEN : BSE1
	J P 8	SHORT: NONE OPEN : S.CHANGE




SEQ PCB

PCB NAME	SW, VR NAME	FUNCTION
HR CRT	V R 1 V R 2 V R 3 V R 4 V R 5	FOCUS
		
	V R 6	Correction of magnification error on CRT center and edge (X scale)
	V R 7	Correction pincushion distortion (X scale)
	V R 8	Correction of magnification error on CRT center and edge (Y scale)
	V R 9	Correction pincushion distortion (Y scale)
	V R 1 0	SPOT POSITION(X scale)
	V R 1 1	SPOT POSITION(Y scale)
	V R 1 2	SQUARING WAVEFORM GAIN (X)
	V R 1 3	SQUARING WAVEFORM GAIN (Y)
	S W 1	SINGLE UNIT ADJUSTMENT




HR PCB

SG Monitor (E)




VIDEO AMP

VIDEO SIGNAL




SH-SBC (1)  
SH-MB (2)




IMEM 1 ②

(1) 特能 TV, 7 回



## PHOTO SIZE



### PHOTO

- PHOTO CRT DEF. AREA :  $126 \times 114\text{mm}$  ( $528 \times 480\text{pixel}$ )
- IMAGE MEMORY AREA :  $122 \times 114\text{mm}$  ( $512 \times 480\text{pixel}$ )
- CRT MONITOR VIEW AREA :  $122 \times 95\text{mm}$  ( $512 \times 400\text{pixel}$ )
- PHOTO FILM AREA :  $114 \times 89\text{mm}$  ( $479 \times 374\text{pixel}$ )
- ADJUSTMENT BOX :  $112 \times 87\text{mm}$  ( $470 \times 365\text{pixel}$ )

### VIEW MONITOR


- CRT MONITOR VIEW AREA :  $195 \times 152\text{mm}$  ( $512 \times 400\text{pixel}$ )
- CRT MONITOR MENU SIZE :  $244 \times 188\text{mm}$  ( $640 \times 480\text{pixel}$ )

VIEW RATIO of MONITOR :  $\times 1.6$

PIXEL ASPECT : 1 : 1

PHOTO PIXEL RESOLUTION : 4.2dot/mm (107DPI)


PCB NAME	SW, VR NAME	FUNCTION
ANALOG	V R 1	SLOW3,4:SG-X CYCLE (Thx=4ms,TP3)
	V R 2	SLOW2 :SG-X CYCLE (Tx=4ms,TP3)
	V R 3	SLOW1 :SG-X CYCLE (Tx=0.7ms,TP3)
	V R 4	SLOW1-4:SG-Y AMPLITUDE(9.1V,TP4)
	V R 5	SPLIT-DUALMAG:SG-X AMPLITUDE
	V R 6	∞ -SCAN:X AMPLITUDE(12V,TP5)
	V R 7	X-RAY MODE:X POSITION
	V R 8	∞ -SCAN:Y AMPLITUDE(8.8V,TP6)
	V R 9	X-RAY MODE:Y POSITION
	V R 1 0	PHOTO-CRT:X WIDTH
	V R 1 1	PHOTO-CRT:X POSITION
	V R 1 2	PHOTO-CRT:Y WIDTH
	V R 1 3	PHOTO-CRT:Y POSITION
	V R 1 4	X-RAY MODE:X AMPLITUDE
	V R 1 5	X-RAY MODE:Y AMPLITUDE
	V R 1 6	(MAG-Y:UPD)
	V R 1 7	MAG-X(COLUMN SCAN X WIDTH)
	V R 1 8	COLUMN SCAN RIGHT ANGLE(ORTHOGONAL)
	V R 1 9	(COLUMN SCAN RIGHT ANGLE:UPD)
	V R 2 0	(MAG-X:UPD)
	V R 2 1	MAG-Y(COLUMN SCAN Y WIDTH)
	V R 2 2	TV-SCAN:X AMPLITUDE
	V R 2 3	TV-SCAN:X POSITION
	V R 2 4	TV-SCAN:Y AMPLITUDE
	V R 2 5	TV-SCAN:Y POSITION
	V R 2 6	ABC-CONTRAST:ADC LEVEL
	V R 2 7	ABC-BRIGHTNESS:ADC LEVEL
	V R 2 8	LINE-ANALYSIS1 AMPLITUDE
	V R 2 9	LINE-ANALYSIS1 POSITION
	V R 3 0	LINE-ANALYSIS2 AMPLITUDE
	V R 3 1	LINE-ANALYSIS2 POSITION
	V R 3 2	TV-SCAN:MEMORY AREA(Thx=10 μS,TP24)
	V R 3 3	PHOTO-CRT BRIGHTNESS(COURSE)
	V R 3 4	PHOTO-CRT BRIGHTNESS(FINE)
	V R 4 0	BSE1 SIGNAL:GAIN
	J P 1	SHORT: NONE OPEN : DBC (X-SCAN)
	J P 2	SHORT: NONE OPEN : DBC (Y-SCAN)
	J P 3	SHORT: NONE OPEN : RASTER ROTATION (X-SCAN)
	J P 4	SHORT: NONE OPEN : RASTER ROTATION (Y-SCAN)
	J P 5	SHORT: NONE OPEN : DBC (BLK)
	J P 6	SHORT: NONE OPEN : DBC (VIDEO)
	J P 7	ALWAYS SHORT




NOTE • at DBC connected, remove jumper pin JP1, JP2, JP6, JP5.

Reserve

PCB NAME	SW, VR NAME	FUNCTION
DEF / LENS	VR 1	OBJ-SPAN
	VR 3	COND1-SPAN
	VR 4	COND2-SPAN
	VR 5	TV-X : MAG CENTER
	VR 6	SLOW-X : MAG CENTER
	VR 7	TV-Y : MAG CENTER
	VR 8	SLOW-Y : MAG CENTER
	VR 9	DYNAMIC FOCUS CENTER
	VR 10	REFERENCE : +10V(TP1)





DEF/LENS PCB




DCPS PCB

PCB NAME	SW,VR NAME	FUNCTION	
S H - S B C	SW 1 A	0 : N-SEM & W 1 : N-SEM & LaBe 2 : H-SEM & W 3 : H-SEM & LaBe 4 : UPD & W 5 : UPD & LaBe 6 : 7 :	8 : 9 : A : B : C : D : E : F :
	SW 1 B	0 : 1 : 2 : 3 : 4 : 5 : 6 : 7 :	8 : 9 : A : B : C : D : E : F :
	SW 2	OFF 1 JAPANESE 2 NTSC 3 NONE 4 NONE 5 RESERVE 6 RESERVE 7 8 NONE	ON ENGLISH PAL X-RAY MODE PSEUDO COLOR LANGUAGE CHANGE
	SW 3	1 NONE 2 NONE 3 NONE 4 RESERVE 5 0-270Pa 6 NONE 7 NONE 8 NONE	RR,DF,TC RS-232C SCSI 0-530Pa BSE2 ELEMENT MOD GRAPHIC
	SW 4	1 RESERVE 2 RESERVE 3 RESERVE 4 5 RESERVE 6 7 SYSTEM 1 OFF 8 SYSTEM 2 OFF	SYSTEM 1 ON SYSTEM 2 ON
	J P 1	SHORT: WDT NORMAL(AUTO RESET) OPEN : WDT NONE (for DEBUG)	
	J P 2	SHORT: BUZZER LOUD OPEN : BUZZER LOW	




PCB NAME	S, VR NAME	FUNCTION
HVC 2	VR 1 VR 2 VR 3	ACC. VOLTAGE ADJUST FILAMENT CURRENT ADJUST FREQUENCY (32 $\mu$ S)




HVC 2 PCB

PCB NAME	SW, VR NAME	FUNCTION
SEQ 2	V R 1	GUN-VACUUM LEVEL ADJ.
	V R 2	CHAMBER-VACUUM LEVEL ADJ.
	V R 3	SPECIMEN-CHANGER VACUUM ADJ.
	J P 1	1-2 SHORT : HIGH 2-3 SHORT : NATURAL
	J P 2	SHORT: NONE OPEN : UTW
	J P 3	1-2 SHORT: W-GUN 2-3 HSORT: LaB6-GUN
	J P 4	1-2 SHORT: W-GUN 2-3 SHORT: LaB6-GUN
	J P 5	ALWAYS SHORT
	J P 6	SHORT: NONE OPEN : BSE2
	J P 7	SHORT: NONE OPEN : BSE1
	J P 8	SHORT: NONE OPEN : S.CHANGE




SEQ PCB

PCB NAME	SW, VR NAME	FUNCTION
HR CRT	VR 1	FOCUS
	VR 2	"
	VR 3	"
	VR 4	"
	VR 5	"
	VR 6	Correction of magnification error on CRT center and edge (X scale)
	VR 7	Correction pincushion distortion (X scale)
	VR 8	Correction of magnification error on CRT center and edge (Y scale)
	VR 9	Correction pincushion distortion (Y scale)
	VR 1 0	SPOT POSITION(X scale)
	VR 1 1	SPOT POSITION(Y scale)
	VR 1 2	SQUARING WAVEFORM GAIN (X)
	VR 1 3	SQUARING WAVEFORM GAIN (Y)
	SW 1	SINGLE UNIT ADJUSTMENT




HR PCB

HS  
'5.8.10  
伊東




鏡体部の内部構造





Angle of Specimen Tilt



Objective Lens of S-2460N




X-ray Analysis of S-3200N



Comparison of Resolution  
S-3200N vs S-2460N




Comparison of Resolution  
S-3200N vs S-2460N




6'1

	Std	Exp 1	Exp 2	Exp 3	備考
A	EDX	CL	CL		
B	エキスレーブ	ジ-ズ	ロビンソン		* デジタル-レギュレーティング・ヒート. N-SHE
C	SED	SED	SED		AUTOMATIC LINE-SUSPEND
D	RXランプ	-	-		
E	(ロック)	(ロック)	(ロック)		
F	(ロック)	(ロック)	(ロック)		
G	EXランプ	EXランプ	EXランプ		* クライオ-ランプ(SP-EX用)
H	ロビンソン	ロビンソン	ロビンソン		UP-TK or SP-EX
I	-	-	CL	-	WDX * I-E-ECI規格品 #フランジ部
J	SP-EX	SP-EX	SP-EX		



8'1(大直径)

	Std	Exp 1	Exp 2	Exp 3	Exp 4	備考
A	CL	EDX	CL	CL	EDX or CL	* BA-1E035規格品. N-SHE
B		ロビンソン	エキスレーブ	ロビンソン	ロビンソン	自動遮断機能
C	SED	SED	SED	SED	SED	
D					(Tニュラ)	* DS-1E7ニンジ. I-E-ECI規格品 (N-SHE)
E	(ロック)	(ロック)	(ロック)	(ロック)	(ロック)	
F	EXランプ	EXランプ	EXランプ	EXランプ	EXランプ	
G	ロビンソン	アニュラ	ロビンソン	ロビンソン	ロビンソン	
H	SP-EX	UP-TK	UP-TK or SP-EX	UP-TK or SP-EX	UP-TK or SP-EX	II-E-1E2ニンジ規格品. #フランジ部
I	-	CL	-	WDX		* I-E-ECI規格品 #フランジ部




A	EDX検出器用ポート
B	ロビンソン反射電子検出器
C	二次電子検出器
D	四分割半導体検出器
E	ステージロック機構
F	試料交換室ランプ
G	空きポート
H	試料交換装置
I	WDX検出器用ポート


大型試料室ポートレイアウト

Objective Lens of S-2460N


30°



X-ray Analysis of S-3200N



Electron Optics of S-3200N



Angle of Specimen Tilt