

# JEM-2100F

*Field Emission Electron Microscope*




# For Nanoanalysis: Field Emi

*Innovative scientific technologies used in today's nanotechnology and bioscience are making remarkable progress. In research on novel materials such as carbon nanotubes, semiconductors and ceramics, as well as basic research in medicine and biology, nanometer-scale inspection, evaluation and analysis are essential.*

*JEOL has developed a new, field emission transmission electron microscope (FE-TEM), the JEM-2100F, a modern instrument for the 21st century.*

*The JEM 2100F integrates a wealth of hardware and software and efficiently acquires valuable research information ... images, diffraction patterns, elemental measurements to name a few ... and enables you to directly view and study nano-scale structures.*



## **Features**

The JEM-2100F is equipped with a field emission electron gun (FEG) that produces high brightness (100 times greater than LaB<sub>6</sub> tip) and is highly stable. This feature is essential for nano-scale ultrahigh resolution and analysis.

The JEM-2100F incorporates multiple additional functions such as a high sensitivity scanning transmission electron microscope (STEM) image observation device and an energy dispersive X-ray spectrometer (EDS) system allowing simple and integrated data acquisition. All operations of the main instrument, including data display and instrument control, are integrated with a PC, thus greatly improving operability.

The EDS system uses a newly developed detector that provides sensitivity more than three times (compared to JEOL conventional detectors) for X-ray low energy regions.

This improvement allows faster analysis with higher reliability.

This combination results in optimum stage movement over the entire magnification range and can be used to stop drift during data acquisition. These benefits lead to efficient imaging and analysis.

# *Transmission Electron Microscope*

## **JEM-2100F**



# Integrated Operation Achieves Excellent Operability

## Multi-functional control

The control system of the JEM-2100F consists of the main system and a number of subsystems. The main system controls basic functions such as the electron optical system and evacuation system. The subsystem controls the electron gun, goniometer stage and panel operation. You can operate the two systems through a host PC built in the instrument designed to comprehensively manage data.


Also, the control system can easily add new functions to the instrument, for example an EDS system and a newly developed scanning image observation device. Furthermore, using the networking capabilities of a PC, you can easily connect various analytical attachments and imaging devices to the instrument.

The control technology of the JEM-2100F integrates a wide range of systems, and is the keystone for adding expanding new TEM application fields.




Analytical capabilities


TEM





STEM & HAADF




EDS



EELS



# Nano Control Mechanism




## New goniometer stage

The newly developed goniometer stage achieves precise specimen movement and stable specimen holding. In particular, the specimen moves with the highest accuracy over a wide range of magnifications: from low magnification to the maximum magnification for nano-area imaging.

This new stage comes with a control mechanism that uses a piezoelectric element, allowing fine specimen movement in nanometer areas.

The stage also incorporates a touch-sensor safety mechanism (JEOL patent), allowing maximum tilt over the eucentric range X, Y and Z movement.




## Large solid angle EDS (option)

To achieve high energy resolution and high sensitivity, the new EDS detector is designed for large solid angle and take-off angle, both of which are essential parameters for X-ray acquisition.

### JED-2300T


Solid angle:	0.28sr
Take-off angle:	24.1°
Detector:	50mm <sup>2</sup>
	At high resolution configuration



## Vibration-proof instrument structure

Measures to prevent vibration are important factors.

The JEM-2100F uses a new base frame and a passive air mount greatly reducing vibration sensitivity.



### Active air mount (option)

The active air mount, having excellent anti-vibration performance, detects floor vibrations and drives an actuator in the opposite direction of the vibrations.



# STEM Expands Analysis Applications



## High resolution STEM images (option)

Control of the scanning image observation device is integrated into the main instrument computer. By the employment of a full digital system, you can easily obtain a variety of high resolution STEM images. The images are displayed on the monitor of the main instrument.


## HAADF\* (option)

The dark-field detector used for imaging is annular in shape. This imaging method reduces the effect of the thickness of the specimen and eliminates diffraction contrast. High-angle incoherent scattered electrons produce images whose contrast depends on Z (atomic number).

\*High-Angle Annular Dark-Field.




## STEM bright filed image (BFI)




HAADF

Specimen: catalyst



HAADF

Specimen: catalyst



Dumbbell image

Specimen: Si (110)

# High Sensitivity High Resolution EDS Analysis

## EDS system (option)





The JED-2300T energy dispersive X-ray spectrometer is a data management system fully integrated into a TEM.

You can obtain the microscope parameters such as magnification and accelerating voltage from the TEM during acquisition.

You can also preset acquisition conditions for elemental mapping and multi-point analysis using STEM image on the screen. In addition, Analysis Station automates acquisition, storage, and management of data.


Together with the STEM drift correction system, Analysis Station enables you to perform high-precision elemental mapping and multi-point analysis because this system allows long-time acquisition with high sensitivity and high energy resolution.

Specimen: carbon nanotube on boron nitride  
Courtesy of Dr. Yoshio Bando  
Advanced Materials Laboratory, National  
Institute for Materials Science.




## Better resolution at low energies

JEOL's unique ultra-thin window detector makes it possible for the JED-2300T to detect light elements. Also, since the JED-2300T offers high energy resolution in low-energy regions, it clearly separates peaks of light elements such as boron (B), carbon (C), nitrogen (N) and oxygen (O). Optionally, the JED-2300T can be equipped with a new detector for high-sensitivity analysis, having a wide effective detection area of 50 mm<sup>2</sup> and an acceptance solid angle of 0.28 sr.



## High energy resolution

The JED-2300T uses a newly developed digital pulse processor and an improved detector. These features allow efficient acquisition of high energy resolution X-ray spectra and high-speed analysis.











# Drift-free Mapping, Spectrum Imaging

## Drift correction system (option)

The STEM drift correction system allows long-time, stable acquisition, making it possible to acquire elemental maps with high spatial resolution and high S/N ratio. In addition, this system provides high energy resolution and high sensitivity to light elements, enabling you to obtain accurate distribution of light elements such as oxygen (O) and nitrogen (N). These features enhance the precision and efficiency of materials evaluation.


- Specimen: DRAM
- Acquisition time: 60 minutes




## EDS active mapping function (spectrum imaging: option)

All pixels in the X-ray mapping data contain spectral information.

After acquiring the mapping data, you can select and add elements to display, set the ROI again, and extract spectra from the desired position (point and area).



Spectra extracted from Area a in the STEM BEI.







Spectra extracted from Area b in the STEM BEI.



# Nano-Area Structure Analysis (EELS Analysis)


Optional EELS maps also allow high sensitivity, high spatial resolution analysis by the use of a high stability, high brightness FEG.

## HAADF images and EELS spectra




Specimen: 64MB DRAM


Spectrum from Point 1



Spectrum from Point 2



Spectrum from Point 3



GATAN EELS spectrometer is installed.

# Automatic 3D Structure Analysis

## TEM 3D tomograph (option)

3D structure analysis macromolecule materials and biological specimens enables one to obtain unique information.

In the TEM observation mode, the specimen is tilted from  $+60^\circ$  to  $-60^\circ$  in fine steps and the recorded images are PC-processed and are subjected to 3D reconstruction for analysis.


### Specimen: Gyroid

Thin section (150 nm to 200 nm)

Tilt angle:  $-60$  to  $+60^\circ$  ( $2.5^\circ$  steps)


Accelerating voltage: 200 kV

Specimen courtesy of Dr. H. Hasegawa of Kyoto University and Dr. K. Jinnai of Kyoto Institute of Technology



## 3D reconstruction system (option)

The specimen is tilted and the images acquired at each step are accumulated. The resulting images are then processed in order to reconstruct a 3D image. This process can be automated on line by linking the high-end algorithm developed by JEOL with a high-precision goniometer.



# High Resolution Supported by High Stability Instrument

## Ultrahigh resolution data showing overall performance


Specimen:  
Amorphous Au particles on Ge film

Direct magnification:  
 $\times 250,000$

Enlarged magnification:  
 $\times 250,000$

Point resolution:  
0.19 nm



Ultrahigh resolution configuration



## Make the spherical aberration coefficient (Cs) = zero: much higher resolution!

You can obtain ultimate high resolution by installing the optional Cs corrector into the JEM-2100F.


\*If you want to use this Cs corrector, please order it together with your JEM-2100F. It cannot be installed at your installation site.




Cs: 0.003 mm  
Point resolution: 0.13 nm

The radius of the circle corresponds to 0.13 nm resolution.


## High resolution image of Si [110]



(Cs=0.03mm)  $\Delta F=-10\text{nm}$



(Cs=-0.03mm)  $\Delta F=10\text{nm}$



Since Cs can be set to minus value, a structure image can be easily obtained even when the specimen is thick and the effects of nonlinear terms cannot be neglected.




# Highly Stable HT System and Electron Optical System

## Highly stable emission current


Highly stable emission of an electron beam, as well as high stability of the high voltage (HT) system, is very crucial to analytical TEM.

Stability of the probe current is essential to achieve reproducibility results during analysis. The figure confirms the stability of the FEG of the JEM-2100F.



## Quick Beam Select System

The Quick Beam Select System is provided for the following four illumination conditions: TEM image observation, EDS analysis, nano-beam diffraction (NBD) and convergent-beam diffraction (CBD). You can select the optimum condition from these four modes with one touch operation. This system enables you to easily perform high precision, microarea analysis, parallel illumination for imaging and CBD.



# Optional Attachments

## Scanning image observation device Full digital control

The BFI and DFI-HAADF detectors for STEM images and the image acquisition unit are provided.

- Scanning image observation device (ASID) EM-24540
- Image acquisition unit EM-27100
- Dark field image observation device EM-24560

## Dry pump system

- TMP rough pump unit EM-23085
- TMP vacuum pump unit EM-23075

## Imaging 3-fold astigmatism correction unit

The components of three-fold astigmatism may have a negative effect on analyzing crystal boundary structures in a high-resolution image or in determining crystal symmetries.

This unit corrects these undesired components and eliminates the effects of three-fold astigmatism.

## Motor drive aperture Remote control

- Motor drive CL aperture EM-20330
- Motor drive OL aperture EM-07320
- Motor drive OL high contrast aperture EM-20340
- Motor drive IL aperture EM-20350
- Motor drive hard X-ray aperture EM-20360
- High tilt retainer EM-21310

## Extensive specimen holders

With the employment of the 5-axis motor-driven microactive goniometer, the tip of the specimen holder is not subject to compressive stress due to atmospheric pressure, thus making it possible for the JEM-2100F to use an ultra-thin specimen holder. The quick specimen retainer, provided as standard, enables you to quickly exchange the specimen only by switching the tip of the holder. Also, extensive specimen holders are available, including the specimen tilting holder, the beryllium specimen tilting holder, the specimen heating holder, and the specimen cooling holder.

## Energy dispersive X-ray spectrometer JED-2300T

Unique selection of detectors with up to 0.28 sr solid angle Excellent resolution at low energies. High speed analysis for sensitive specimen characterization Drift free mapping for light element mapping with high spatial resolution.





EU's product safety standard



Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

# JEOL

JEOL Ltd. 1-2 Musashino 3-chome Akishima Tokyo 196-8558 Japan Sales Division ☎(042)528-3381 〓(042)528-3386

<http://www.jeol.co.jp/>

## ARGENTINA

COASIN S. A. C. I. yF.  
Virrey del Pino 4071,  
1430 Buenos Aires, Argentina  
Telephone: 54-11-4552-3185  
Facsimile: 54-11-4555-3321

## AUSTRALIA & NEW ZEALAND

JEOL (AUSTRALASIA) Pty. Ltd.  
Unit 9, 750-752 Pittwater Road,  
Brookvale, N. S. W. 2100, Australia  
Telephone: 43-2233-53838  
Facsimile: 61-2-9905-8255

## AUSTRIA

LABCO GmbH  
Dr.-Trittemmel-Gasse 8,  
A-3013 Pressbaum, Austria  
Telephone: 43-2233-53838  
Facsimile: 61-2-9905-8286

## BANGLADESH

A.Q. CHOWDHURY & CO. Pvt. Ltd.  
Baridhara Central Plaza 87, Suhrawardy Avenue  
2nd Floor Baridhara,  
Dhaka-12129 Bangladesh  
Telephone: 880-2-9862272, 9894583  
Facsimile: 880-2-988070

## BELGIUM

JEOL (EUROPE) B. V.  
Zaventem/Ikaros Business Park,  
Ikaroslaan 7A, B-1930 Zaventem, Belgium  
Telephone: 32-2-720-0560  
Facsimile: 32-2-720-0134

## BRAZIL

FUGIWARA ENTERPRISES  
INSTRUMENTOS CIENTIFICOS LTDA.  
Avenida Iteberaba,3563  
02739-000 Sao Paulo, SPI Brazil  
Telephone: 55-11-3983-8144  
Facsimile: 55-11-3983-8140

## CANADA

JEOL CANADA, INC.  
(Represented by Soquelec, Ltd.)  
5757 Cavendish Boulevard, Suite 540,  
Montreal, Quebec H4W 2W8, Canada  
Telephone: 1-514-482-6427  
Facsimile: 1-514-482-1929

## CHILE

TECSIS LTDA.  
Avenida Holanda 1248,  
Casilla 50/9 Correo 9, Providencia, Santiago, Chile  
Telephone: 56-2-205-1313  
Facsimile: 56-2-225-0759

## CHINA

JEOL LTD., BEIJING OFFICE  
Room No. B2101/2106,  
Vantone New World Plaza,  
No. 2 Fuwai Street, Xicheng District,  
Beijing 100037, P. R. China  
Telephone: 86-10-6804-6321/6322/6323  
Facsimile: 86-10-6804-6324

## JEOL LTD., SHANGHAI OFFICE

Sanhe Building 11F2, Yan Ping Road,  
No. 121, Shanghai 200042, P.R. China  
Telephone: 86-21-6246-2353  
Facsimile: 86-21-6246-2836

## JEOL LTD., GUANG ZHOU OFFICE

(South Power) World Trade Center Building  
371-375, Huang Shi East-Road, Guang Zhou,  
510095, P. R. China  
Telephone: 86-20-8778-7848  
Facsimile: 86-20-8778-4268

## CYPRUS

MESLO LTD.  
Scientific & Laboratory Division,  
P. O. Box 27709, Nicosia Cyprus  
Telephone: 357-2-666070  
Facsimile: 357-2-660355

## EGYPT

JEOL SERVICE BUREAU  
3rd Fl. Nile Center Bldg.,  
Nawal Street, Dokki, (Cairo), Egypt  
Telephone: 20-2-335-7220  
Facsimile: 20-2-338-4186

## FRANCE

JEOL (EUROPE) S. A.  
Espace Claude Monet,  
1, Allée de Giverny 78290  
Croissy-sur-Seine, France  
Telephone: 33-13015-3737  
Facsimile: 33-13015-3747

## GERMANY

JEOL(GERMANY)GmbH  
Oskar-Von-Miller-Strasse 1,  
85386 Eching Germany  
Telephone: 49-8165-77346  
Facsimile: 49-8165-77512

## GREAT BRITAIN & IRELAND

JEOL (U.K.) LTD.  
JEOL House, Silver court, Watchmead, Welwyn  
Garden City, Herts AL7 1LT., England  
Telephone: 44-1707-377117  
Facsimile: 44-1707-373254

## GREECE

N. ASTERIAS S. A.  
56-58, S. Trikoupi Str. P.O.Box 26140  
GR-10022 Athens, Greece  
Telephone: 30-1-823-5383  
Facsimile: 30-1-823-9567

## HONG KONG

FARMING LTD.  
Unit 1009, 10/F., MLC Millennium Plaza  
663 King's Road, North Point, Hong Kong  
Telephone: 852-2815-7299  
Facsimile: 852-2581-4635

## INDIA

Blue Star LTD. (HQ)  
Analytical Instruments Department,  
'Sahas'414/2 Veer Savarkar Marg,  
Prabhadery Mumbai 400 025, India  
Telephone: 91-22-5666-4068  
Facsimile: 91-22-5666-4001

## Blue Star LTD. (New Delhi)

Analytical Instruments Department,  
E-44/12 Okhla Industrial Area,  
Phase-11, New Delhi 110 020, India  
Telephone: 91-11-5149-4000  
Facsimile: 91-11-5149-4004

## Blue Star LTD. (Calcutta)

Analytical Instruments Department,  
7, Hare Street Calcutta 700 001, India  
Telephone: 91-33-2248-0131  
Facsimile: 91-33-2248-1599

## Blue Star LTD. (Chennai)

Analytical Instruments Department, Lakshmi  
Neela Rite Choice Chambers, 5, Bazullah Road,  
3rd Floor T. Nagar Chennai 600 017, India  
Telephone: 91-44-2815-8846  
Facsimile: 91-44-2815-8015

## INDONESIA

PT. TEKNOLABindo PENTA PERKASA  
J1. Gading BukitRaya,  
Komplek Gading Bukit Indah Blok I/11,  
Kelapa Gading Jakarta 14240, Indonesia  
Telephone: 62-21-45847057/58/59  
Facsimile: 62-21-45842729

## IRAN

IMACO LTD.  
No. 141 Felestin Ave.,  
P. O. Box 13145-537, Tehran, Iran  
Telephone: 98-21-6402191/6404148  
Facsimile: 98-21-8978164

## ITALY

JEOL (ITALIA) S.p.A.  
Centro Direzionale Green Office Via Dei Tulpiani,  
1, 20090 Pieve, Emanuele (MI), Italy  
Telephone: 39-2-9041431  
Facsimile: 39-2-90414353

## KOREA

JEOL KOREA LTD.  
Sunmin Bldg. 6th F1.,218-16, Nonhyun-Dong,  
Kangnam-Ku, Seoul, 135-010, Korea  
Telephone: 82-2-511-5501  
Facsimile: 82-2-511-2635

## KUWAIT

YUSUF I. AL-GHANIM & CO. (YIACO)  
P. O. Box 435, 13005 - Safat, Kuwait  
Telephone: 965-4832600/4814358  
Facsimile: 965-4844954/4833612

## MALAYSIA

JEOL (MALAYSIA) SDN. BHD. (359011-M)  
205, Block A, Mezzanine Floor,  
Kelana Business Center97,  
Jalan SS 7/2, Kelana Jaya,  
47301 Petaling Jaya, Selangor, Malaysia  
Telephone: 60-3-7492-7722  
Facsimile: 60-3-7492-7723

## MEXICO

JEOL DE MEXICO S.A. DE C.V.  
Av. Amsterdam #46 DEPS. 402  
Col. Hipodromo, 06100 Mexico D.F. Mexico  
Telephone: 52-5-55-211-4511  
Facsimile: 52-5-55-211-0720

## PAKISTAN

ANALYTICAL MEASURING SYSTEM (PVT) LTD.  
(AMS LTD.)  
F-8/1-5, KDA Scheme No. 1,  
Main Tipu Sultan Road,  
P. O. Box 10604, Karachi-75350, Pakistan  
Telephone: 92-21-4525185/4525186  
Facsimile: 92-21-4316331/4316332

## PANAMA

PROMED S.A.  
Parque Industrial Costa del Este  
Urbanizacion Costa del Este  
Apartado 6281, Panama, Panama  
Telephone: 507-269-0044  
Facsimile: 507-263-5622

## PHILIPPINES

PHILAB INDUSTRIES INC.  
7487 Bagtikan Street, SAV Makati, 1203 Metro,  
Manila Philippines  
Telephone: 63-2-896-7218  
Facsimile: 63-2-897-7732

## PORTUGAL

Izasa. Portugal Lda.  
R. do Proletariado 1,  
2790-138 CARNAXIDE Portugal  
Telephone: 351-21-424-7300  
Facsimile: 351-21-418-6020

## SAUDI ARABIA

ABDULREHMAN ALGOSAIBI G. T.B.  
Algosaihi Bldg., Airport Rd., P. O. Box 215,  
Riyadh 11411, Saudi Arabia  
Telephone: 966-1-479-3000  
Facsimile: 966-1-477-1374

## SCANDINAVIA

JEOL (SKANDINAVISKA) A.B.  
Hammarbacken 6 A, Box 716  
191 27 Sollentuna, Sweden  
Telephone: 46-8-28-2800  
Facsimile: 46-8-29-1647

## SERVICE & INFORMATION OFFICE

JEOL NORWAY  
Ole Deviks vei 28, N-0614 Oslo, Norway  
Telephone: 47-2-2-64-7930  
Facsimile: 47-2-2-65-0619

## JEOL FINLAND

Ylakaupinkuja 2, FIN-02360 Espoo, Finland  
Telephone: 358-9-8129-0350  
Facsimile: 358-9-8129-0351

## JEOL DENMARK

Naverland 2, DK-2600 Glostrup, Denmark  
Telephone: 45-4345-3434  
Facsimile: 45-4345-3433

## SINGAPORE

JEOL ASIA PTE. LTD.  
29 International Business Park,  
#04-02A Acer Building,  
Tower B Singapore 609923  
Telephone: 65-6565-9989  
Facsimile: 65-6565-7552

## SOUTH AFRICA

ADI Scientific (Pty) Ltd.  
109 Blandford Road, North Riding,Randburg  
(PO box 71295 Bryanston 2021)  
Republic of South Africa  
Telephone: 27-11-462-1363  
Facsimile: 27-11-462-1466

## SPAIN

IZASA. S.A.  
Aragoneses, 13,  
28100 Alcobendas,  
(Poligono Industrial) Madrid, Spain  
Telephone: 34-91-663-0500  
Facsimile: 34-91-663-0545

## SWITZERLAND

JEOL(GERMANY)GmbH  
Oskar-Von-Miller Strasse 1,  
85386 Eching Germany  
Telephone: 49-8165-77346  
Facsimile: 49-8165-77512

## TAIWAN

JIE DONG CO., LTD.  
7F, 112, Chung Hsiao East Road, Section 1, Taipei,  
Taiwan 10023, Republic of China  
Telephone: 886-2-2395-2978  
Facsimile: 886-2-2322-4655

## JEOL TAIWAN SEMICONDUCTORS LTD.

11F, No. 346, Pei-Ta Road, Hsin-Chu City 300,  
Taiwan Republic of China  
Telephone: 886-3-523-8490  
Facsimile: 886-2-523-8503

## THAILAND

BECTHAI BANGKOK EQUIPMENT  
& CHEMICAL CO., Ltd.  
300 Phaholyothin Rd. Phayathai,  
Bangkok 10400, Thailand  
Telephone: 66-2-615-2929  
Facsimile: 66-2-615-2350/2351

## THE NETHERLANDS

JEOL (EUROPE) B.V.  
Tupolevlaan 28-A, 1119 NZ Schiphol-Rijk,  
The Netherlands  
Telephone: 31-20-6533088  
Facsimile: 31-20-6531328

## TURKEY

TEKSER KOLL. STI.  
Acibadem Cad. Erdem Sok, Bayer Art, N° 6/1,  
81010, Uskudar, Istanbul, Turkey  
Telephone: 90-216-3274041  
Facsimile: 90-216-3274046

## UAE

BUSINESS COMMUNICATIONS LLC.  
P. O. Box 2534, Abu Dhabi UAE  
Telephone: 971-2-6348495  
Facsimile: 971-2-6316465

## USA

JEOL USA, INC.  
11 Dearborn Road, Peabody, MA. 01960, U. S. A.  
Telephone: 1-978-535-5900  
Facsimile: 1-978-536-2205/2206

## JEOL USA, INC. WEST OFFICE

5653 Stoneridge Drive Suite  
#110 Pleasanton, CA. 94588 U. S. A.  
Tel: 1-925-737-1740  
Fax: 1-925-737-1749

## VENEZUELA

MITSUBISHI VENEZOLANA C. A.  
Avenida Francisco de Miranda Los Palos Grandes,  
Caracas, Venezuela  
Telephone: 58-212-209-7402  
Facsimile: 58-212-209-7496